

MRI Scan-Time Reduction Using Complex Valued AI In Image Space

by

Md Tauheed Ansari
(21AG3AI06)

Under the guidance of
Prof. Mahesh Mohan M R




Department Artificial Intelligence
Indian Institute of Technology, Kharagpur

Contents




1.	Introduction
2.	Justification
3.	Objectives
4.	Literature Review
	- MRI Fundamentals
	- Existing Methods
5.	Materials And Methods
6.	Results and Discussions
7.	Conclusion
8.	Future works and References



INTRODUCTION

- Magnetic Resonance Imaging (MRI)
 - medical imaging technique
 - visualizes anatomy and physiological processes
 - without ionizing radiation.
- MRI uses strong magnetic fields, gradients, and radio waves to generate detailed images.
- Unlike CT and PET scans, MRI does not involve harmful radiation, making it safer for repeated diagnostics.
- MRI offers superior soft tissue contrast, useful in
 - brain
 - muscle
 - abdominal imaging.




JUSTIFICATION

- **Long scan times** cause discomfort for injured and paining patients.
- **Motion artifacts** reduce image quality during lengthy scans.
- **Higher operational costs** due to low patient throughput.
- **Limited accessibility** in busy or resource-constrained hospitals.
- **Need for faster MRI** to improve patient care and system efficiency.





OBJECTIVES

- **Reduce MRI scan time** without compromising image quality.
 - **Enhance diagnostic throughput** by accelerating MRI imaging.
 - **Apply complex-valued AI techniques** for faster k-space data reconstruction.
 - **Evaluate performance** using standard image quality metrics (e.g., PSNR, SSIM).
 - **Develop a scalable solution** for real-time MRI scan acceleration.
- 

LITERATURE REVIEW

Research Paper	Year	Major Findings
SENSE Pruessmann et al.	1999	Introduced sensitivity encoding (SENSE) , enabling faster MRI scans by exploiting multiple receiver coils for parallel imaging, thus reducing scan times.
GRAPPA Griswold et al.	2002	Developed GRAPPA (Generalized Autocalibrating Partially Parallel Acquisition), a method that accelerates MRI by using parallel coils and calibration-free reconstruction.
ESPIRiT Uecker et al.	2014	Proposed ESPIRiT (Extended Sensitivity Encoding with Iterative Reconstruction), enhancing parallel MRI techniques by improving coil sensitivity maps for more accurate reconstructions.
Compressed Sensing Jaspan et el.	2015	Applied compressed sensing to MRI, showing how under-sampled data can still lead to high-quality images, significantly reducing scan time without sacrificing image resolution.
Deep SENSE Xi Peng et al.	2021	Combined deep learning with SENSE to improve image reconstruction quality and accelerate MRI scans, providing a hybrid approach that outperforms traditional methods.

MRI Fundamentals

MRI Steps



1. Magnetization

- Strong magnetic field aligns body's H protons
- Protons align either parallel or antiparallel to the field, creating a net magnetization.



2. Radiofrequency (RF) Pulses

- RF pulse is applied to disturb proton alignment
- This energy tips protons into a higher energy state.



3. Signal Detection

- Once the RF pulse is switched off, protons 'relax'
- They emit radio signals that are detected by coils around body.

MRI Fundamentals

MRI Steps





4. K-Space Sampling




5. Image Reconstruction

- The signals are collected in K-space (a frequency domain).
- Each point in K-space holds phase and frequency information needed to reconstruct the final image.
- Using mathematical algorithms (like Fourier Transform), the data from K-space is processed into a detailed 2D or 3D image.




K-Space Visualisations



Aliasing Problem

Aliasing in Accelerated MRI


- Acceleration skips k-space lines \rightarrow undersampling.
- Causes aliasing: overlapping image structures.
- Multicoil data helps resolve aliasing using the concept of coil sensitivity.
- Reconstruction methods (e.g., SENSE) correct the artifacts.



Parallel Imaging - Our best tool

Two kinds of Techniques-



1. **IMAGE- DOMAIN** based technique
-performed in image space *after* reconstruction of data from the individual coils
2. **K-SPACE** based technique
-operate primarily on *k*-space data *before* image reconstruction




SENSE (1999)

Sensitivity Encoding for Fast MRI


Visual representation of Sampling under Acceleration




(a) Random mask with 4-fold acceleration



(b) Random mask with 8-fold acceleration






(c) Equispaced mask with 4-fold acceleration



(d) Equispaced mask with 8-fold acceleration

ALIASING- Resolution Approach





Materials and Methods

- **fastMRI Dataset: An Open Dataset and Benchmarks for Accelerated MRI**

Designed to advance MRI reconstruction.

The fastMRI raw k-space data is stored in the HDF5 (.h5) file format, which efficiently organizes large, complex datasets, enabling fast access to MRI raw data and associated metadata for training and evaluation.

Data Types:

Raw multi-coil k-space data (complex-valued measurements)

Emulated single-coil k-space data (simulated from multi-coil)

Ground-truth images (root-sum-of-squares reconstruction)

DICOM images (diverse scanners, post-processed)

Tasks Enabled:

Single-coil reconstruction

Multi-coil reconstruction

Modalities & Organs:

Knee MRI (1,594 multi-coil raw scans; 10,000 DICOM scans)

Brain MRI (6,970 multi-coil raw scans; 10,000 DICOM scans)

BART: Berkeley Advanced Reconstruction Toolbox

-Computational Magnetic Resonance Imaging

- Free and open-source framework for **Computational MRI**.
- Designed **for research use only** (not for diagnostic use).
- Over **100 commands** covering diverse reconstruction tasks.
- Strong support for **FFT** and **Non-Uniform FFT (NUFFT)** operations.
- Facilitates **parallel imaging** (SENSE, POCSENSE) and **compressed sensing**.
- Scriptable via **command-line interface** for easy automation.

The logo for BART (Berkeley Advanced Reconstruction Toolbox) features the word "BART" in a large, bold, black, sans-serif font. The letters are enclosed within a white rectangular frame that has a stylized, slightly irregular border, resembling a thick L-shaped bracket or a frame with a small gap at the top-left and bottom-right corners.

Key Features-

- Nonlinear inversion and compressed sensing reconstruction.
- Deep learning integration with TensorFlow and Reconet.
- Synthetic data generation (phantoms, fake k-space).
- Image evaluation tools (NMRSE, ROI stats, MIP views).
- Wavelet transforms, TGV and ICTV regularization.
- Fast, lightweight, suitable for rapid MRI research prototyping.

BART. Available commands are:

avg	bench	bin	bitmask	cabs	calc
caldir	calmat	carg	casorati	cc	ccapply
cdf97	circshift	conj	conv	conway	copy
cpyphs	creal	crop	delta	ecalib	ecaltwo
epg	estdelay	estdims	estshift	estvar	extract
fakeksp	fft	fftmod	ffttrot	fftshift	filter
flatten	flip	fmac	fovshift	homodyne	ictv
index	invert	itsense	join	looklocker	lrmatrix
mandelbrot	measure	mip	mnist	moba	mobafit
morphop	multicfl	nlinv	nlmeans	nnet	noise
normalize	nmse	nufft	nufftbase	onehotenc	ones
pattern	phantom	pics	pocsense	poisson	pol2mask
poly	psf	reconet	repmat	reshape	resize
rmfreq	rof	roistat	rss	rtnlinv	sake
saxpy	scale	sdot	show	signal	sim
slice	spow	sqpics	squeeze	ssa	std
svd	tensorflow	tg	threshold	toimg	traj
transpose	twixread	upat	var	vec	version
walsh	wave	wavelet	wavepsf	whiten	window
wshfl	zeros	zexp			

Methodology: MRI Reconstruction Pipeline

1. Input Handling

- Supports both `.cfl` (BART-native) and `.h5` (fastMRI) data formats.
- Data is preprocessed and formatted for compatibility with BART tools.

2. Preprocessing

- `.cfl`: Loaded directly with `readcfl`, already compatible with BART.
- `.h5`: Slice selected, transposed, reshaped, and written to `.cfl` using `writetcfl`.

3. Coil Sensitivity Map Estimation

- Used BART's `ecalib` command for ESPIRiT calibration.
- `-r24` selects central k-space lines; `-m1` generates one map per coil for SENSE.

4. Undersampling via Masking

- **fastMRI-style Variable-Density Masks**: Created using `RandomMaskFunc`, controlled by `center_fraction` and `acceleration`.
- **Custom Centered Masks**: Manually defined masks for clearer control and interpretability

5. Image Reconstruction

- **Reference:** Inverse FFT followed by coil combination using **fft** and **fmac**.
- **SENSE:** BART's **pics** with regularized least-squares and single-set sensitivity maps.
- **ESPIRiT:** Robust multi-map method using **pics** with ESPIRiT sensitivity maps.

6. Visual Evaluation

- Side-by-side comparison of reference, SENSE, and ESPIRiT reconstructions.
- Focused on aliasing, blurring, and sharpness.

7. Quantitative Analysis

- Generated **difference maps** to measure deviation from the reference.
- Highlighted areas of high error, especially at tissue boundaries or low-SNR regions.

DATALOADER

Dataset Preparation

- Defined `Custom_FastMRIDataset` class based on `torch.utils.data.Dataset`.
- Loaded `.npy` MRI data with `mmap_mode='r'` to save RAM.
- Applied undersampling masks on k-space (`apply_mask`).

DataLoader Setup

- Created `train_loader_check` with batch size 4 and shuffling.
- Checked data shapes and types to verify correct loading.

Transformations

- Converted k-space to image-space (`convert_K_to_I`).
- Created Root-Sum-of-Squares (RSS) combined images

Flexible Data Return

- Return options controlled by flags: k-space (`K`), image-space (`I`), or RSS-combined (`rss_combine`)

Visualization

- Picked first item from batch (`batchidx=0`).
- Displayed the ground truth (`loader_rss`) and masked version (`loader_masked_rss`) using `matplotlib.pyplot.imshow()` with grayscale colormap.

Results and Discussions

BART Implementation

- Implemented BART on CFL files
- Implemented BART on fastMRI data file

1. Reference Scan Visualization

Downsampled
Scan with AF=2




Calibration: `bart ecalib`
estimates coil sensitivity maps from
k-space data.

Fourier Transform: `bart fft`
`-i` converts k-space data to coil
images.

Combination: `bart fmac -C`
performs a **sensitivity-weighted
sum** across coils to generate the
final image.

Scan Reconstruction Result Comparison

1. Reference Scan



2. ESPIRiT



3. SENSE




4. Downsampled(with AF=2)



Difference Map Comparison

1. Reference Scan




2. ESPIRiT



3. SENSE




4. Downsampled(with AF=2)




Scan Reconstruction Result Comparison


1. Reference Scan




2. ESPIRiT



3. SENSE




4. Downsampled(with AF=4)



Difference Map Comparison


1. Reference Scan



2. ESPIRiT



3. SENSE




4. Downsampled(with AF=4)




ESPIRiT Results for Various Acceleration Factors


Acceleration Factor=2



Acceleration Factor=3




Acceleration Factor=4



Thus, a need for better Deep Learning alternatives arises.

CVNN Training Loss Result Plot



- Sharp drop in **training loss** during early epochs.
- **Validation loss** decreases steadily, indicating good generalization.
- Both losses **stabilize after ~5 epochs**, showing convergence.
- No overfitting observed – validation closely follows training loss.
- Confirms effective training of the complex-valued U-Net for 2x upscaling.

Conclusion

Implemented a complete MRI reconstruction pipeline:

- Used *BART toolbox* to reconstruct MR images with classical methods like **SENSE**, **ESPIRiT**, and **PICS**.
- Understood and applied concepts of fourier transforms, *coil sensitivity estimation*, and *parallel imaging*.

Built robust masking and preprocessing functions:

1. Designed custom **undersampling masks** including acceleration factor and centre fraction as arguments, suitable for accelerated MRI acquisition simulation.
2. Developed *efficient preprocessing steps* like **k-space masking**, **inverse Fourier transforms**, **root sum-of-squares (RSS) image combination**.
3. Carefully structured the pipeline to maintain consistency between masked and unmasked data pairs.

Designed robust masking functions and preprocessing pipelines (k-space masking, Fourier transforms, RSS combination) ensuring clean and standardized inputs.

Trained and Plotted training and validation loss curves for complex-valued **U-Net** to verify effective learning and convergence during 2x upscaling training.

Built flexible, memory-efficient DataLoaders supporting multiple data formats (k-space, image-space, RSS) for scalable training.

Future Works

1. **Extend current classical reconstruction methods** by incorporating *Compressed Sensing MRI* techniques for better sparsity-based recovery.
2. **Implement deep learning-driven reconstruction models** like *DeepSense* and *SENSE with GAN* for learned k-space to image-space mappings.
3. **Explore advanced architectures** such as *Self-Attention CNNs* and *High-Frequency Diffusion Models* to enhance structural preservation in reconstructions.
4. **Develop densely connected and transformer-based models** (*DenseNet*, *Swin Transformer*, *ReconFormer*) for faster and more accurate MRI reconstruction.
5. **Benchmark and compare all DL models** against the strong *mathematical baselines* (*SENSE*, *ESPIRiT*, *PICS*) already established.

References used:

- ❑ Hirose, A. *Complex-Valued Neural Networks: The Merits and Their Origins*, IJCNN, IEEE, 2009.
- ❑ Dash, B.N., Khare, N. *Deep Complex Neural Network Applications in Remote Sensing*, SPIE, 2021.
- ❑ Dola, A.D. *Complex Valued CNNs for Image Restoration*, MSc Thesis, IIT Kharagpur, 2023.
- ❑ Zbontar, J. et al. *fastMRI: An Open Dataset and Benchmarks for Accelerated MRI*, arXiv:1811.08839, 2018.
- ❑ Pruessmann, K.P. et al. *SENSE: Sensitivity Encoding for Fast MRI*, MRM, 1999.
- ❑ Uecker, M. et al. *ESPIRiT: An Eigenvalue Approach to Parallel MRI*, MRM, 2014.
- ❑ Peng, X. et al. *DeepSENSE: Learning Coil Sensitivity with Deep Learning*, MRM, 2022.
- ❑ Lv, J. et al. *GANs for Multi-Channel MRI Reconstruction*, Comp. in Biology & Medicine, 2021.
- ❑ Lee, C. et al. *Complex-Valued Neural Networks: A Comprehensive Survey*, IEEE/CAA JAS, 2022.



Thanks

&

Questions are Welcomed!

